



QUICKIE[®]
Live without limits.

Argon[®]

Fine tune your
performance...

Quickie Argon[®]

Finely and easily adjustable

"Fine tune" your performance ... first time, second time, every time! Infinite fine tuning possibilities to ensure perfect driving ergonomics can be achieved and maintained.



Fine tune your performance



Excellent handling!

High tolerances and rigidity for excellent handling, rolling performance and cornering. Rigid frame construction and excellent driving characteristics. Minimalist open frame design ensures easy lifting and transportation.



True open frame design



Rear axle with infinite rear wheel centre of gravity adjustment for optimum driving and turning performance



Lightweight axle plate and axle tube clamp offers independent adjustments of seat height and axle spacing



Improved, durable folding backrest system offering angle adjustment from -12° to $+8^{\circ}$

Looks to perform...

and a custom choice to hit the road. A wide range of options, six different upholstery colours and 20 brilliant and matt finishes for your frame, hubs, forks, wheels, handrims and footplates.



Lightweight tubular footrest (option)



Compact wheel lock folds away nicely below the frame (option)



Fold down push handles (option)




Fold down backrest with "double lock" (option) holds the backrest down for easy transportation

Argon



* Based on 0° camber and 24" wheels



Technical data sheet

20 individual frame colours	Weight of heaviest part	7kg
6 individual sling colours	Turning circle	max 115cm (45.2")
Total Weight	at least 9.9kg	Transit tested to ISO 7176 Part 19 and passed
Maximum user weight	120kg (18.8 stone)	Yes

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.

Find out more today!

Visit www.sunrisemedical.com for more information, to find your local dealer or to request a demonstration.



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Argon and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairsUK



Sunrise Medical Limited

Thorns Road,
Brierley Hill
West Midlands
DY5 2LD

Tel: +44 (0) 845 605 66 88

Fax: +44 (0) 845 605 66 89

www.SunriseMedical.com