# Quickie®

**Ultra-Light** 

The Quickie 2 is already the most popular ultra-light folding frame wheelchair available today, and we've found a way to make it even better. With twelve new features and options, the Quickie 2 now offers more fine tune adjustments for improved versatility, customization and ease of use. With the Quickie 2 Family the possibilities are endless.



#### Quickie 2 Highlights:

- Fully adjustable axle plate with camber, center of gravity, wheel base width and seat height adjustments
- Even more caster adjustability to accommodate any modifications of the frame that affect alignment
- Angle-adjustable backrest
- Fully modular design for complete customization
- The Quickie True Fit Program allows the chair to grow
- High quality Quickie components and hardware ensure durability and performance
- Widest variety of hanger, footplate, frame, caster, wheel, backrest and armrest options
- New Jay Precision™ Back option
- Travel wheel option
- One-arm drive option (see Quickie Kids picture)
- Uni-lateral (one-hand) wheel lock system option







# LDING FAMILY

### Quickie Folding Chairs – A complete folding line to suit your mobility needs.

### **Ultra-Light**













Lightweight









Positioning

**Heavy Duty** 







## Quickie 2 & Quickie 2 Kids & Quickie 2 HD

Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	a rall	5" LP POLY 7" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PN 5 <sup>1/4"</sup> FORK
Тор	16 <sup>15/16</sup> "	18 <sup>3/4</sup> "	18 <sup>5/16</sup> "	18 <sup>3/16</sup> "	18 <sup>1/2</sup> "
Bottom	17 15/16"	19 3/4"	19 <sup>5/16</sup> "	19 3/16"	19 <sup>1/2</sup> "

#### Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only) 20" 22" 22" 24" 24" 24" 24" 26" AXLE/HOLE PNEU FP PNEU LP FP HP PNEU/

POSITION		POLY		POLY	POLY	CLNCHR	KEVLAR	CLNCHR	
1 AND 2	18 <sup>7/8</sup> "	19 <sup>1/2</sup> "	19 3/4"	20"	20 1/4"	20 3/8"	20 1/2"	20 1/4"	
2 AND 3	17 3/4"	18 <sup>3/8"</sup>	18 <sup>5/8</sup> "	18 <sup>7/8"</sup>	19 <sup>1/8"</sup>	19 <sup>1/4</sup> "	19 <sup>3/8"</sup>	20 1/8"	
3 AND 4	16 <sup>5/8</sup> "	17 <sup>1/4</sup> "	17 <sup>1/2</sup> "	17 <sup>3/4</sup> "	18"	18 <sup>1/8"</sup>	18 <sup>1/4</sup> "	19"	
4 AND 5	15 1/2"	16 1/8"	16 3/8"	16 5/8"	16 7/8"	17"	17 1/8"	17 7/8"	

#### Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH 12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
SPOKE WHEEL 201/2"	211/2"	221/2"	231/2"	241/2"	251/2"	261/2"	271/2"	281/2"	291/2"	301/2"
W/ HANDRIMS										

Quickie 2	Quickie 2 HD	Quickie 2 Kids

Weight Capacity: 250 lbs Weight Capacity: 350 lbs Weight Capacity: 250 lbs

#### Quickie 2HP

Front Seat Height (For Standard Frames Only)

CASTER FORK	3"MICROS	3"MICROS	3"MICROS	5"LP POLY	5"LP POLY	5"LP PC
MOUNTING HOLE	5 <sup>1/4</sup> "FORK	6"FORK	7"FORK	5 <sup>1/4</sup> "FORK	6"FORK	7"FORK
TOP	16 <sup>5/16</sup> "	17 1/8"	18 <sup>1/8</sup> "	17 3/16"	18"	19"
воттом	17 5/16"	18 <sup>1/8</sup> "	19 <sup>1/8</sup> "	18 3/16"	19"	20"

#### Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR	22"	24"	24"	24"	24" HT	26"	26" HT
AXLE/HOLE	PNEU	LATEX	HP	PNEU	TEN/B-BALL	LATEX	PNEU
POSITION		TUBULAR	CLNCHR		TURBO TRD	HP CLNCHR	
1 AND 2	19 3/4"	20 1/4"	20 3/8"	20 1/2"	20 5/8"	21 1/4"	21 5/8"
2 AND 3	18 <sup>5/8</sup> "	19 <sup>1/8</sup> "	19 <sup>1/4</sup> "	19 <sup>3/8</sup> "	19 <sup>1/2</sup> "	20 1/8"	20 1/2"
3 AND 4	17 1/2"	18"	18 <sup>1/8</sup> "	18 <sup>1/4</sup> "	18 <sup>3/8</sup> "	19"	19 3/8"
4 AND 5	16 <sup>3/8</sup> "	16 <sup>7/8</sup> "	17"	17 <sup>1/8</sup> "	17 1/4"	17 1/8"	18 <sup>1/4</sup> "

#### Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	1/"	18"	19"	20"	
SPOKE/MAG WHEEL W/HANDRIMS	21 1/2"	22 1/2"	231/2"	241/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"	

#### Quickie 2 HP

Weight: 22 lbs Weight Capacity: 250 lbs

#### Quickie M6

FRAME WIDTH:	24"	26"	28"	30"
OVERALL WIDTH:	32.5"	34.5"	36.5"	38.5
FOLDED WIDTH:	14"			

#### Quickie M6

Weight: 53 lbs. w/legrests (26 x 20)

Weight Capacity: 650 lbs

Qu	ickie	Rec	liner
_			

ront	Seat	Heig	h
------	------	------	---

	9					
CASTER FORK	5" LP POLY	6" PNEU	6" POLY	8" PNEU	8" POLY	
MOUNTING HOLE	6" FORK	6" FORK	6" FORK	5 <sup>1/4</sup> " FORK	5 <sup>1/4"</sup> FORK	
TOP	17 3/4"	18 <sup>5/16</sup> "	18 <sup>3/16</sup> "	N/A	N/A	
воттом	18"	19 <sup>5/16</sup> "	19 3/16"	19 <sup>1/2</sup> "	19 5/16"	

#### **Rear Seat Height**

I riour Cour rioigi					
REAR	22"	24"	24"	24"	
AXLE HOLE	FP	FP	HP	Pneu/	
POSITION	Poly	Poly	Clincher	Kevlar	
Тор	15"	15 <sup>3/4"</sup>	15 <sup>7/8</sup> "	16"	
Middle	16"	16 <sup>3/4"</sup>	16 <sup>7/8"</sup>	17"	
Bottom	17"	17 3/4"	17 7/8"	18"	

#### Armrest Height

Height-Adjustable

9" - 14"

#### Footrest Height

SWING-AWAY	ELR	
14" - 21"	16 1/2" - 19"	

#### **Overall Width (Without Armrests)**

Frame Width	14"	15"	16"	17"	18"	19"	20"	Folded
SPOKE/MAG WHEEL								
W/ HANDRIMS &								
QUICK-RELEASE AXLES	22 1/4"	23 1/4"	24 1/4"	25 <sup>1/4</sup> "	26 <sup>1/4</sup> "	27 1/4"	28 1/4"	12 <sup>1/2</sup> "
SPOKE/MAG WHEEL								
W/ HANDRIMS &								
THREADED AYLES	22"	23"	24"	25"	26"	27"	28"	12 1/2"

Quickie Recliner
Weight: 42 lbs. (w/ swing-away footrests, height-adjustable armrests, headrest extension and mag. wheels

Weight Capacity: 250 lbs

#### Quickie LX & LXI Front Seat Height

CASTER FORK	5"LP Poly	5"LP Poly	6" Poly	6"Pneu	8"Poly	8" Pneu	8"x2"Caster
MOUNTING HOLE	4 <sup>3/4</sup> " Fork	51/4" Fork	6"Fork	6"Fork	51/4"Fork	51/4" Fork	8"x 2"Fork
TOP	N/A	N/A	17"	17"	N/A	N/A	N/A
воттом	16"	16 1/2"	18"	18"	18"	18"	N/A

#### Rear Seat Height

REAR AXLE HOLE	20" LP	20"	22" LP	22"	24" LP	24"	
POSITION (LX ONLY)	POLY	PNEU	POLY	PNEU	POLY	PNEU	
TOP	15"	15 <sup>1/2</sup> "	16"	16 <sup>1/2</sup> "	17"	17 <sup>1/2</sup> "	
воттом	16"	16 <sup>1/2</sup> "	17"	17 <sup>1/2</sup> "	18"	18 <sup>1/2</sup> "	

REAR AXLE HOLE	20" LP	20"	22" LP	22"	24" LP	24"	
POSITION (ADJ.AXLE)	POLY	PNEU	POLY	PNEU	POLY	PNEU	
TOP	N/A	15"	15 <sup>1/2</sup> "	16"	16"	17"	
MIDDLE	15 <sup>1/2</sup> "	16"	16 <sup>1/2</sup> "	171/2"	17 <sup>1/2</sup> "	18"	
воттом	17"	171/2"	18"	181/2"	181/2"	191/2"	

#### Armrest Height

#### **Footrest Height**

FLIP-BACK	FLIP-BACK/	PADDED S/A	HT. ADJ.	SWING-AWAY	ELR/ALR
	HEIGHT ADJ.				
9 1/4"	8 <sup>3/4</sup> "-12 <sup>5/8</sup> "	8 3/4"	7 <sup>1/4</sup> "-12 <sup>1/4</sup> "	12 <sup>3/4</sup> "-15 <sup>1/4</sup> "	14"-16 <sup>1 /2"</sup>
			8 3/4"-13 3/4"	15 <sup>1/4</sup> "-17 <sup>3/4</sup> "	16 <sup>1/2</sup> "-19"
				17 <sup>3/4</sup> "-20 <sup>1/4</sup> "	

#### Overall Width (Without Armrests)

FRAME WIDTH	12"	14"	16"	18"	20"	FOLDED	
MAG WHEEL							
W/ HANDRIMS	19"	21"	23"	25"	27"	121/2"	
SPOKEWHEEL W/ HANDRIMS	20"	22"	24"	26"	28"	121/2"	
W HANDRING	20		<u>_</u> +	20	20	12 -	

Quickie Recliner	
Weight:	28 lbs
Weight Capacity:	250 lbs

# **SPECIFICATIONS**

CASTER FORK	/ICRO	5"	POLY	OLY 6" PNEU						
MOUNTING HOLE		1/4"	5	1/4"	5 1/4"		5 1/4"			
TOP		V/A	17	7 1/4"	17 <sup>1/8</sup>		17 3/4"			
воттом	17	7 1/4"	18	3 1/4"	18 <sup>7/8</sup>		18 3/4"			
REAR SEAT HEIGHT (FOR STANDARD FRAMES AND AXLE POSITION ONLY)										
AXLE HOLE	PNEU	PNEU	LATEX	TEN/BBALL	HP	PNEU	TEN/BBALL			
POSITION	. = 4/41	4 (01)		TURBO TRD			TURBO TRI			
0°,4°/NO SPACER	15 1/4"	16 <sup>1/8</sup> "	15 3/4"	16 1/4"	16 7/8"	17 1/4"	17 5/8"			
8°, 12°/NO SPACER	14 1/2"	15 <sup>3/8</sup> "	15"	15 <sup>1/2</sup> "	16 <sup>1/8</sup> "	16 <sup>1/2</sup> "	16 <sup>7/8</sup> "			
8°, 4° + 3/4°	16"	16 <sup>7/8</sup> "	16 <sup>1/2</sup> "	17"	17 <sup>5/8</sup> "	18"	18 <sup>3/8</sup> "			
8°, 12° + 3/4°	15 <sup>1/4</sup> "	16 <sup>1/8</sup> "	15 <sup>3/4</sup> "	16 <sup>1/4</sup> "	16 <sup>7/8</sup> "	17 <sup>1/4</sup> "	17 <sup>5/8</sup> "			
0°, 4° + 1 1/2°	16 <sup>3/4</sup> "	17 <sup>5/8</sup> "	17 1/4"	17 3/4"	18 3/8"	18 <sup>3/4</sup> "	19 <sup>1/8</sup> "			
8°, 12° + 1 1/2°	15"	16 <sup>7/8</sup> "	16 <sup>1/2</sup> "	17"	17 <sup>5/8</sup> "	18"	18 <sup>3/8</sup> "			
Overall Width	•			11						
Frame Width	14"	15"	16"	17"	18"	19"	20"			
SpokeWheel*										
with Handrims	24"	25"	26"	27"	28"	29"	30"			
Quickie M6										
Weight:		23 lbs.								

#### Quickie 2 True Fit Program

Physical changes can cause a chair not to fit properly. To ensure proper fit, Sunrise offers one-time replacement for growth of the seat sling, back upholstery, footplates, hangers, back tubes, frame and crossbrace at no extra charge. (For up to three years from delivery on all Quickie 2 family wheelchairs shipped after November 15, 1995.)

ACTUAL MEASUREMENTS MAY VARY ACCORDING TO CHAIR CONFIGURATION AND ACCESSORIES. PLEASE CONSULT ORDER FORM FOR COMBINATIONS AND RESTRICTIONS. SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE.

LIFETIME FRAME AND AXLE WARRANTY

We highly recommend distribution through a certified Rehab Professional when fitting and ordering any Quickie Products.

